

# ADAM8 antibody



Catalog Number: 107768

## Product name

ADAM8 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ADAM8 Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells, mouse pancreas tissue. Positive IHC detected in human tonsillitis tissue, human pancreas tissue. Observed molecular weight by Western-blot: 65 kDa

## Preparation

This antibody was obtained by immunization of ADAM8 recombinant protein (Accession Number: NM\_001109). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:107768(ADAM8 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107768(ADAM8 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog

# ADAM8 antibody



Catalog Number: 107768

---

No:107768(ADAM8 Antibody) at dilution of 1:50 (under 40x lens)