

# ADCY1 antibody

Catalog Number: 107779

---

## Product name

ADCY1 antibody

## Specificity

Human; other species not tested.

## Antibody description

ADCY1 Rabbit Polyclonal antibody. Positive WB detected in Y79 cells. Observed molecular weight by Western-blot: 123kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_021116). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:107779(ADCY1 antibody) at dilution of 1:500