

ADCY1 antibody



Catalog Number: 107779

Product name

Polyclonal

ADCY1 antibody

Ig Type

Rabbit IgG

Specificity

Human; other species not tested.

Applications

ELISA, WB

Antibody description

ADCY1 Rabbit Polyclonal antibody. Positive WB detected in Y79 cells. Observed molecular weight by Western-blot: 123kd

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_021116).
Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:107779(ADCY1 antibody) at dilution of 1:500