

ABCF1 antibody



Catalog Number: 107805

Product name

ABCF1 antibody

IP: 1:200-1:2000

IF: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ABCF1 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive IP detected in K-562 cells. Positive WB detected in K-562 cells. Observed molecular weight by Western-blot: 110 kDa

Preparation

This antibody was obtained by immunization of ABCF1 recombinant protein (Accession Number: NM_001025091). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IP, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:107805(ABCF1 antibody) at dilution of 1:1500



IP Result of anti-ABCF1 (IP:Catalog No:107805, 4ug; Detection:Catalog No:107805 1:500) with K-562 cells lysate 1200ug.



Immunofluorescent analysis of (-20oc Ethanol) fixed HeLa cells using Catalog No:107805(ABCF1

ABCF1 antibody



Catalog Number: 107805

Antibody) at dilution of 1:50 and Alexa Fluor 488-
conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)