

ABCG1 antibody

Catalog Number: 107807

Product name

ABCG1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ABCG1 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in multi-cells/tissue, HeLa cells, Jurkat cells, mouse brain tissue, mouse lung tissue, RAW264.7 cells. Positive IHC detected in human cervical cancer tissue, human tonsillitis tissue. Observed molecular weight by Western-blot: 60-65kd

Preparation

This antibody was obtained by immunization of ABCG1 recombinant protein (Accession Number: NM_016818). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


IP: 1:500-1:5000

IHC: 1:20-1:200

Validations




WB result of Catalog No:(ABCG1 antibody) on multiple lysates.




Immunohistochemical of paraffin-embedded human cervical cancer using Catalog No:107807(ABCG1 antibody) at dilution of 1:100 (under 10x lens)

ABCG1 antibody

Catalog Number: 107807



Immunohistochemical of paraffin-embedded human cervical cancer using Catalog No:107807(ABCG1 antibody) at dilution of 1:100 (under 40x lens) 

IP Result of anti-ABCG1 (IP:Catalog No:107807, 4ug; Detection:Catalog No:107807 1:1000) with mouse brain tissue lysate 4000ug.