

ABHD15 antibody



Catalog Number: 107813

Product name

ABHD15 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ABHD15 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 52 kDa

Preparation

This antibody was obtained by immunization of ABHD15 recombinant protein (Accession Number: NM_198147). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:107813(ABHD15 Antibody) at dilution of 1:300