

ABHD15 antibody



Catalog Number: 107813

Product name

Polyclonal

ABHD15 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse; other species not tested.

Applications

Antibody description

ELISA, WB

ABHD15 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 52 kDa

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of ABHD15 recombinant protein (Accession Number: NM_198147). Purification method: Antigen affinity purified.

Formulation


PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:107813(ABHD15 Antibody) at dilution of 1:300