

# ABT1 antibody



Catalog Number: 107821

## Product name

ABT1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ABT1 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HeLa cells. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 31-32 kDa

## Preparation

This antibody was obtained by immunization of ABT1 recombinant protein (Accession Number: NM\_013375). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107821(ABT1 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human brain using Catalog No:107821(ABT1 antibody) at dilution of 1:100 (under 10x lens)