

ABT1 antibody



Catalog Number: 107821

Product name

ABT1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ABT1 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HeLa cells. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 31-32 kDa

Preparation

This antibody was obtained by immunization of ABT1 recombinant protein (Accession Number: NM_013375). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107821(ABT1 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human brain using Catalog No:107821(ABT1 antibody) at dilution of 1:100 (under 10x lens)