

# ACAA2 antibody

Catalog Number: 107824

## Product name

ACAA2 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

ACAA2 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer. Positive WB detected in HeLa cells, HEK-293 cells. Observed molecular weight by Western-blot: 42kd

## Preparation

This antibody was obtained by immunization of ACAA2 recombinant protein (Accession Number: NM\_006111). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107824(ACAA2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:107824(ACAA2 antibody) at dilution of 1:100 (under 10x lens)