

# AGFG1 antibody

Catalog Number: 107826

## Product name

AGFG1 antibody

WB: 1:200-1:2000

IP: 1:200-1:1000

## Specificity


Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

## Antibody description

AGFG1 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, HeLa cells, HepG2 cells, HT-1080 cells, Jurkat cells, K-562 cells. Positive IP detected in mouse testis tissue. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 58 kDa

## Validations




## Preparation

This antibody was obtained by immunization of AGFG1 recombinant protein (Accession Number: NM\_004504). Purification method: Antigen affinity purified.

mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107826(AGFG1 antibody) at dilution of 1:400

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:107826(AGFG1 antibody) at dilution of 1:100 (under 10x lens)

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications

ELISA, WB, IHC, IP

## Dilutions

Recommended Dilution:



# AGFG1 antibody



Catalog Number: 107826

---

IP Result of anti-AGFG1 (IP:Catalog No:107826,

4ug; Detection:Catalog No:107826 1:400) with mouse testis tissue lysate 4000ug.