

AGGF1 antibody



Catalog Number: 107828

Product name

AGGF1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

AGGF1 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HEK-293 cells, HeLa cells, rat liver tissue. Positive IHC detected in human spleen tissue. Observed molecular weight by Western-blot: 84-100 kDa

Preparation

This antibody was obtained by immunization of AGGF1 recombinant protein (Accession Number: BC029382). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:10-1:100

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107828(AGGF1 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human spleen using Catalog No:107828(AGGF1 antibody) at dilution of 1:25 (under 40x lens)