

# AGPAT1 antibody

Catalog Number: 107831

## Product name

AGPAT1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

AGPAT1 Rabbit Polyclonal antibody. Positive IP detected in Raji cells, mouse skeletal muscle tissue. Positive IHC detected in human skeletal muscle tissue, human testis tissue. Observed molecular weight by Western-blot: 32 kDa

## Preparation

This antibody was obtained by immunization of AGPAT1 recombinant protein (Accession Number: NM\_006411). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, IP

## Dilutions

Recommended Dilution:


IP: 1:200-1:1000

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human skeletal muscle using Catalog No:107831(AGPAT1 antibody) at dilution of 1:50 (under 10x lens)



IP Result of anti-AGPAT1 (IP:Catalog No:107831, 4ug; Detection:Catalog No:107831 1:300) with Raji cells lysate 2000ug.