

**Product name**

AHCY antibody

**Specificity**

Human, Mouse, Rat; other species not tested.

**Antibody description**

AHCY Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue, HeLa cells. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 48 kDa

**Preparation**

This antibody was obtained by immunization of AHCY recombinant protein (Accession Number: NM\_001362750). Purification method: Antigen affinity purified.

**Formulation**

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**

Polyclonal

**Ig Type**

Rabbit IgG

**Applications**


ELISA, WB, IHC

**Dilutions**


Recommended Dilution:

WB: 1:500-1:5000

IHC: 1:20-1:200

**Validations**

mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:107842(AHCY antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:107842(AHCY antibody) at dilution of 1:100 (under 10x lens)