

# AHSA2 antibody



Catalog Number: 107848

## Product name

Polyclonal

AHSA2 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse; other species not tested.

## Applications

## Antibody description

ELISA, WB

AHSA2 Rabbit Polyclonal antibody. Positive WB detected in fetal human brain tissue. Observed molecular weight by Western-blot: 34 kDa

## Dilutions


Recommended Dilution:

WB: 1:200-1:1000

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_152392). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107848(AHSA2 Antibody) at dilution of 1:300