

# AID antibody



Catalog Number: 107849

---

## Product name

Polyclonal

AID antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

AID Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 48-52 kDa

## Dilutions

Recommended Dilution:

## Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of AID recombinant protein (Accession Number: NM\_020661). Purification method: Antigen affinity purified.

## Validations

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107849(AID antibody) at dilution of 1:800

## Clonality