

AIP/ARA9 antibody



Catalog Number: 107856

Product name

AIP/ARA9 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

AIP/ARA9 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, HeLa cells, NIH/3T3 cells. Positive FC detected in HeLa cells. Positive IHC detected in human liver cancer tissue, human heart tissue. Observed molecular weight by Western-blot: 38 kDa

Preparation

This antibody was obtained by immunization of AIP/ARA9 recombinant protein (Accession Number: NM_003977). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, FC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:107856(AIP antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:107856(AIP antibody) at dilution of 1:50 (under 40x lens)




Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:107856(AIP antibody) at dilution of 1:50 (under 10x lens)

AIP/ARA9 antibody



Catalog Number: 107856



1X10⁶ HeLa cells were stained with 0.2ug AIP antibody (Catalog No:107856, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.