

ALKBH4 antibody



Catalog Number: 107857

Product name

ALKBH4 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

ALKBH4 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, K-562 cells, MCF7 cells, mouse ovary tissue. Positive IF detected in HepG2. Observed molecular weight by Western-blot: 34kd

Preparation

This antibody was obtained by immunization of ALKBH4 recombinant protein (Accession Number: NM_017621). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:107857(ALKBH4 antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using ALKBH4 antibody Catalog No:107857 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).