

ALLC antibody

Catalog Number: 107860

Product name

ALLC antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

Validations

Antibody description

ALLC Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, human heart tissue, human skeletal muscle tissue. Positive IHC detected in human testis tissue, human kidney tissue, human spleen tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 46kd



Preparation

This antibody was obtained by immunization of ALLC recombinant protein (Accession Number: BC029652). Purification method: Antigen affinity purified.

mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:107860(ALLC antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


ELISA, WB, IF, IHC

Dilutions

Recommended Dilution:




Immunofluorescent analysis of HepG2 cells, using ALLC antibody Catalog No:107860 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



ALLC antibody

Catalog Number: 107860

Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:107860(ALLC Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:107860(ALLC Antibody) at dilution of 1:50 (under 40x lens)