

alpha Actin antibody

Catalog Number: 107863

Product name

alpha Actin antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

alpha Actin Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue, mouse skeletal muscle tissue, mouse thymus tissue. Positive IP detected in mouse heart tissue. Observed molecular weight by Western-blot: 42 kDa

Preparation

This antibody was obtained by immunization of alpha Actin recombinant protein (Accession Number: NM_001406483). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:107863(ACTC1 antibody) at dilution of 1:800



IP Result of anti-ACTC1 (IP:Catalog No:107863, 4ug; Detection:Catalog No:107863 1:800) with mouse heart tissue lysate 3200ug.