

alpha Actinin antibody



Catalog Number: 107865

Product name

alpha Actinin antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

alpha Actinin Rabbit Polyclonal antibody. Positive WB detected in C6 cells, human spleen tissue, mouse spleen tissue, NIH/3T3 cells. Observed molecular weight by Western-blot: 100-105kd

Preparation

This antibody was obtained by immunization of alpha Actinin recombinant protein (Accession Number: NM_001102). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



C6 cells were subjected to SDS PAGE followed by western blot with Catalog No:107865(a-Actinin Antibody) at dilution of 1:2000