

ALPPL2 antibody



Catalog Number: 107882

Product name

ALPPL2 antibody

ELISA, IHC

Dilutions

Recommended Dilution:

Specificity


Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

Antibody description

ALPPL2 Rabbit Polyclonal antibody. Positive IHC detected in human placenta tissue, human small intestine tissue.

Validations




Preparation

This antibody was obtained by immunization of ALPPL2 recombinant protein (Accession Number: NM_031313). Purification method: Antigen affinity purified.

Immunohistochemistry of paraffin-embedded human placenta tissue slide using Catalog No:107882 (ALPPL2 Antibody) at dilution of 1:200 (under 10x lens).

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Immunohistochemistry of paraffin-embedded human placenta tissue slide using Catalog No:107882 (ALPPL2 Antibody) at dilution of 1:200 (under 40x lens).

Applications