

AF-6 antibody



Catalog Number: 107900

Product name

Polyclonal

AF-6 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

AF-6 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, A431 cells, human brain tissue. Observed molecular weight by Western-blot: 182-210kd

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_005936). Purification method: Antigen affinity purified.

Validations



Formulation

0.1M NaHCO₃, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:107900(MLLT4 antibody) at dilution of 1:500