

# AFAP1 antibody

Catalog Number: 107902

## Product name

AFAP1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

AFAP1 Rabbit Polyclonal antibody. Positive IF detected in MCF-7 cells. Positive WB detected in rat brain tissue, DU 145 cells, HeLa cells, MCF7 cells, mouse brain tissue. Observed molecular weight by Western-blot: 110kd

## Preparation

This antibody was obtained by immunization of AFAP1 recombinant protein (Accession Number: NM\_001371090). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107902(AFAP1 Antibody) at dilution of 1:1000



Immunofluorescent analysis of MCF-7 cells using Catalog No:107902(AFAP1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)