

AFF1-Specific antibody



Catalog Number: 107905

Product name

Polyclonal

AFF1-Specific antibody

Ig Type

Rabbit IgG

Specificity

human; other species not tested.

Applications

ELISA, WB

Antibody description

AFF1-Specific Rabbit Polyclonal antibody. Positive WB detected in K-562 cells. Observed molecular weight by Western-blot: 125-145 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_005935). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:107905(AFF1-Specific antibody) at dilution of 1:300