

AFF4 antibody



Catalog Number: 107907

Product name

AFF4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

AFF4 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, HeLa cells, HepG2 cells. Positive IP detected in HepG2 cells. Observed molecular weight by Western-blot: 140kd

Preparation

This antibody was obtained by immunization of AFF4 recombinant protein (Accession Number: BC063007). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:107907(AFF4 antibody) at dilution of 1:500



IP Result of anti-AFF4 (IP:Catalog No:107907, 4ug; Detection:Catalog No:107907 1:300) with HepG2 cells lysate 3200ug.