

AIRE antibody



Catalog Number: 107920

Product name

AIRE antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

AIRE Rabbit Polyclonal antibody. Positive IHC detected in human stomach tissue. Positive WB detected in L02 cells, human plasma tissue, human spleen tissue, mouse thymus tissue, U-937 cells. Observed molecular weight by Western-blot: 65-70 kDa

Preparation

This antibody was obtained by immunization of AIRE recombinant protein (Accession Number: NM_000658). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:107920(AIRE Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human stomach slide using Catalog No:107920(AIRE Antibody) at dilution of 1:50