

AKR1D1 antibody



Catalog Number: 107944

Product name

AKR1D1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

AKR1D1 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 33 kDa

Preparation

This antibody was obtained by immunization of AKR1D1 recombinant protein (Accession Number: NM_005989). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107944(AKR1D1 Antibody) at dilution of 1:600