

# ALDH1L1 antibody



Catalog Number: 107964

## Product name

ALDH1L1 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ALDH1L1 Rabbit Polyclonal antibody. Positive IP detected in mouse liver tissue. Positive WB detected in mouse liver tissue. Positive IHC detected in mouse brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 100 kDa

## Preparation

This antibody was obtained by immunization of ALDH1L1 recombinant protein (Accession Number: NM\_144776). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:200-1:1000

## Validations




mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107964(ALDH1L1 Antibody) at dilution of



1:300

IP Result of anti-ALDH1L1 (IP:Catalog No:107964, 4ug; Detection:Catalog No:107964 1:500) with mouse liver tissue lysate 4000ug.



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog

# ALDH1L1 antibody



Catalog Number: 107964

---

No:107964(ALDH1L1 Antibody) at dilution of 1:50 (under 40x lens)