

ALDH8A1 antibody

Catalog Number: 107975

Product name

ALDH8A1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

ALDH8A1 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, mouse kidney tissue. Positive IHC detected in human hepat°Cirrhosis. Observed molecular weight by Western-blot: 53kd

Preparation

This antibody was obtained by immunization of ALDH8A1 recombinant protein (Accession Number: NM_022568). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:107975(ALDH8A1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:107975(ALDH8A1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:107975(ALDH8A1 antibody) at dilution of 1:50 (under 40x lens)