

ALY antibody



Catalog Number: 107986

Product name

Polyclonal

ALY antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

ALY Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells. Observed molecular weight by Western-blot: 32 kDa

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:1000

This antibody was obtained by immunization of ALY recombinant protein (Accession Number: BC052302). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107986(ALY Antibody) at dilution of 1:300