

AMIGO3 antibody



Catalog Number: 107994

Product name

AMIGO3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

AMIGO3 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, HeLa cells, human brain tissue, L02 cells, mouse kidney tissue, mouse ovary tissue. Observed molecular weight by Western-blot: 45-50 kDa

Preparation

This antibody was obtained by immunization of AMIGO3 recombinant protein (Accession Number: NM_198722). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107994(AMIGO3 antibody) at dilution of 1:300