

# AMOT antibody

Catalog Number: 107998

## Product name

AMOT antibody

## Specificity

Human; other species not tested.

## Antibody description

AMOT Rabbit Polyclonal antibody. Positive IF detected in HEK-293 cells. Positive IP detected in HEK-293 cells. Positive WB detected in HEK-293. Observed molecular weight by Western-blot: 80 kDa; 130 kDa

## Preparation

This antibody was obtained by immunization of AMOT recombinant protein (Accession Number: NM\_133265). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IP, WB, IF


## Dilutions

Recommended Dilution:

WB: N/AIP: 1:200-1:1000


IF: 1:20-1:200

## Validations




WB: anti-Amot 24550 antibody 1:1000

WB results of Catalog No:107998 (1:1000) with over-expressed AMOT by Customer.



IP Result of anti-AMOT (IP:Catalog No:107998, 4µg; Detection:Catalog No:107998 1:300) with HEK-293 cells lysate 2500µg.



Immunofluorescent analysis of (-20oc Ethanol) fixed HEK-293 cells using Catalog No:107998(AMOT Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat

# AMOT antibody



Catalog Number: 107998

---

Anti-Rabbit IgG(H+L)