

AMPD1-Specific antibody



Catalog Number: 108001

Product name

Polyclonal

AMPD1-Specific antibody

Ig Type

Rabbit IgG

Specificity

human, mouse, rat; other species not tested.

Applications

ELISA, WB

Antibody description

AMPD1-Specific Rabbit Polyclonal antibody.
Positive WB detected in mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 75-83kd

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000036).
Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:108001(AMPD1-Specific antibody) at dilution of 1:300