

AMPD2 antibody



Catalog Number: 108002

Product name

Polyclonal

AMPD2 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

AMPD2 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue. Observed molecular weight by Western-blot: 101kd

Dilutions

Recommended Dilution:


Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of AMPD2 recombinant protein (Accession Number: NM_139156). Purification method: Antigen affinity purified.

Validations

Formulation



PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:108002(AMPD2 antibody) at dilution of 1:500