

Amylase, alpha antibody



Catalog Number: 108014

Product name

Amylase, alpha antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Amylase, alpha Rabbit Polyclonal antibody.
Positive IHC detected in human pancreas tissue.
Positive WB detected in human liver tissue, mouse pancreas tissue. Observed molecular weight by Western-blot: 58 kDa

Preparation

This antibody was obtained by immunization of Amylase, alpha recombinant protein (Accession Number: NM_000699). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:108014(AMY2A antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human pancreas using Catalog No:108014(AMY2A antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human pancreas using Catalog No:108014(AMY2A antibody) at dilution of 1:100 (under 40x lens)