

APOC4 antibody

Catalog Number: 108029

Product name

APOC4 antibody

Specificity

Human; other species not tested.

Antibody description

APOC4 Rabbit Polyclonal antibody. Positive WB detected in human plasma tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 17 kDa

Preparation

This antibody was obtained by immunization of AP^oC4 recombinant protein (Accession Number: NM_001646). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

Validations



human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:108029(APoc4 Antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells using Catalog No:108029(APoc4 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)