

# APOF antibody

Catalog Number: 108031

## Product name

APOF antibody

## Specificity

human; other species not tested.

## Antibody description

APOF Rabbit Polyclonal antibody. Positive IHC detected in human liver. Positive WB detected in human blood tissue. Observed molecular weight by Western-blot: 29 kDa

## Preparation

This antibody was obtained by immunization of APOF recombinant protein (Accession Number: NM\_001638). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



human blood tissue were subjected to SDS PAGE followed by western blot with Catalog No:108031(APOF antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human liver using Catalog No:108031(APOF antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:108031(APOF antibody) at dilution of 1:50 (under 40x lens)