

APOM antibody

Catalog Number: 108039

Product name

APOM antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

APOM Rabbit Polyclonal antibody. Positive IHC detected in human hepat°Cirrhosis tissue. Positive WB detected in human plasma tissue. Observed molecular weight by Western-blot: 25 kDa

Preparation

This antibody was obtained by immunization of APOM recombinant protein (Accession Number: NM_019101). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:108039(APOM Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human hepatocirrhosis tissue slide using Catalog No:108039(APOM Antibody) at dilution of 1:50 (under 40x lens)