

APPL2 antibody



Catalog Number: 108043

Product name

Polyclonal

APPL2 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

APPL2 Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells. Observed molecular weight by Western-blot: 74 kDa

Dilutions

Recommended Dilution:


Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of APPL2 recombinant protein (Accession Number: NM_018171). Purification method: Antigen affinity purified.

Validations

Formulation



PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:108043(APPL2 Antibody) at dilution of 1:600

Clonality