

androgen receptor antibody



Catalog Number: 108049

Product name

androgen receptor antibody

Specificity

Human, Monkey; other species not tested.

Antibody description

androgen receptor Rabbit Polyclonal antibody.
Positive IP detected in HepG2 cells. Positive WB
detected in HEK-293 cells, COS-7 cells, HepG2
cells. Observed molecular weight by Western-blot:
110 kDa

Preparation

This antibody was obtained by immunization of
androgen receptor recombinant protein
(Accession Number: NM_000044). Purification
method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE
followed by western blot with Catalog
No:108049(AR antibody) at dilution of 1:1000



IP Result of anti-AR (IP:Catalog No:108049, 4ug;
Detection:Catalog No:108049 1:500) with HepG2
cells lysate 3600ug.