

Angiotensin 2 antibody



Catalog Number: 108052

Product name

Angiotensin 2 antibody

WB: 1:200-1:2000

IHC: 1:20-1:200

Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

Antibody description

Angiotensin 2 Rabbit Polyclonal antibody. Positive WB detected in HUVEC cells, K-562 cells, mouse brain tissue. Positive IHC detected in human liver cancer tissue, human breast cancer tissue. Positive IF detected in A549 cells, PC-12 cells. Observed molecular weight by Western-blot: 68 kDa

Validations




Preparation

This antibody was obtained by immunization of Angiotensin 2 recombinant protein (Accession Number: NM_001147). Purification method: Antigen affinity purified.

HUVEC cells were subjected to SDS PAGE followed by western blot with Catalog No:108052(ANGPT2 Antibody) at dilution of 1:600

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:108052(ANGPT2 Antibody) at dilution of 1:50 (under 10x lens)

Ig Type


Rabbit IgG

Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:




Angiopoietin 2 antibody



Catalog Number: 108052

Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:108052(ANGPT2 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of A549 cells using Catalog No:108052(ANGPT2 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)