

ANGPTL7 antibody



Catalog Number: 108060

Product name

ANGPTL7 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ANGPTL7 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, COLO 320 cells, HepG2 cells, Y79 cells. Observed molecular weight by Western-blot: 45 kDa

Preparation

This antibody was obtained by immunization of ANGPTL7 recombinant protein (Accession Number: NM_021146). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:108060(ANGPTL7 antibody) at dilution of 1:800