

ANKRD29 antibody

Catalog Number: 108074

Product name

ANKRD29 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ANKRD29 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 29kDa

Preparation

This antibody was obtained by immunization of ANKRD29 recombinant protein (Accession Number: NM_173505). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108074(ANKRD29 Antibody) at dilution of



1:300

IP Result of anti-ANKRD29 (IP:Catalog No:108074, 3ug; Detection:Catalog No:108074 1:300) with mouse brain tissue lysate 2640ug.