

ANXA2R antibody



Catalog Number: 108110

Product name

ANXA2R antibody

Specificity

Human; other species not tested.

Antibody description

ANXA2R Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer tissue, human brain tissue.

Preparation

This antibody was obtained by immunization of ANXA2R recombinant protein (Accession Number: NM_001382352). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using Catalog No:108110(C5orf39 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using Catalog No:108110(C5orf39 Antibody) at dilution of 1:50 (under 40x lens)