

ANXA9 antibody

Catalog Number: 108111

Product name

ANXA9 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ANXA9 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human liver tissue, human heart tissue, human kidney tissue, human placenta tissue, human skin tissue, human spleen tissue, human testis tissue. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 38 kDa

Preparation

This antibody was obtained by immunization of ANXA9 recombinant protein (Accession Number: BC005830). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:1000


IHC: 1:20-1:200

IF: 1:10-1:100

Validations




mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:108111(ANXA9 antibody) at dilution of 1:300




Immunohistochemical of paraffin-embedded human liver using Catalog No:108111(ANXA9 antibody) at dilution of 1:50 (under 10x lens)

ANXA9 antibody

Catalog Number: 108111



Immunohistochemical of paraffin-embedded human liver using Catalog No:108111(ANXA9 antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using ANXA9 antibody Catalog No:108111 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).