

# AP3B1 antibody

Catalog Number: 108121

## Product name

AP3B1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

AP3B1 Rabbit Polyclonal antibody. Positive IP detected in COLO 320 cells. Positive WB detected in mouse thymus tissue, COLO 320 cells. Observed molecular weight by Western-blot: 140 kDa

## Preparation

This antibody was obtained by immunization of AP3B1 recombinant protein (Accession Number: NM\_003664). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:108121(AP3B1 antibody) at dilution of 1:500



IP Result of anti-AP3B1 (IP:Catalog No:108121, 3ug; Detection:Catalog No:108121 1:500) with COLO 320 cells lysate 2500ug.