

APC1 antibody

Catalog Number: 108132

Product name

APC1 antibody

WB: 1:1000-1:10000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

APC1 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue, human cervical cancer tissue, human stomach tissue. Positive WB detected in HeLa cells, HEK-293 cells, HT-1080 cells, K-562 cells, mouse brain tissue, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 200-210 kDa

Preparation

This antibody was obtained by immunization of APC1 recombinant protein (Accession Number: NM_022662). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108132(ANAPC1 antibody) at dilution of 1:1500



Immunohistochemical of paraffin-embedded human testis using Catalog No:108132(ANAPC1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human testis using Catalog No:108132(ANAPC1

APC1 antibody



Catalog Number: 108132

antibody) at dilution of 1:50 (under 40x lens)