

# APH1A antibody

Catalog Number: 108139

## Product name

APH1A antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

APH1A Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, human brain tissue. Observed molecular weight by Western-blot: 30 kDa

## Preparation

This antibody was obtained by immunization of APH1A recombinant protein (Accession Number: NM\_001077628). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108139(APH1A antibody) at dilution of 1:500



IP Result of anti-APH1A (IP:Catalog No:108139, 4ug; Detection:Catalog No:108139 1:300) with HeLa cells lysate 3200ug.