

APLF antibody



Catalog Number: 108143

Product name

APLF antibody

Specificity

human; other species not tested.

Antibody description

APLF Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells. Observed molecular weight by Western-blot: 60kd

Preparation

This antibody was obtained by immunization of APLF recombinant protein (Accession Number: NM_173545). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:108143(APLF antibody) at dilution of 1:600